

# soluble TNF receptor 1 (sTNF RI) serum

Test Code: 30144

## Clinical and Procedure

### Clinical Utility

For the quantitative measurement of soluble Tumor Necrosis Factor receptor 1 (sTNFR1). The pro-inflammatory cytokine, TNF $\alpha$  and its soluble receptor, sTNFR1, are potent modulators of the inflammatory process.

### Procedure

The assay for quantification of sTNFR1 is a sandwich ELISA performed in a microtiter plate format. Conversion of a chromogenic substrate produces a color, the intensity of which is proportional to the concentration of sTNFR1 in the sample material. A standard curve is used to calculate the concentration of sTNFR1 in each of the test samples. This test has not been cleared or approved for diagnostic use by the U.S. Food and Drug Administration.

### Specificity

Specific to human sTNFR1.

## Turnaround Time

3 business days from receipt of specimen

## Specimen Information

Specimen Type	Order Code	CPT Code	NY Approved	Volume	Assay Range
serum	30144	84238	Yes	1 mL	0.2 - 16 ng/mL

#### Special Instructions

- Whole blood should be collected in serum tube.
- Allow to clot for 30 to 60 minutes and centrifuged to isolate the serum.
- 1 mL of plasma sample should be removed to a sterile tube and frozen immediately (-70°C).

The result is reported in ng/mL. The assay range is approximately 0.2 to 16 ng/mL. The reference range for a healthy population is less than 2.4 ng/mL. However it should be noted that these ranges are obtained from a limited population of apparently healthy adults and are not diagnostic thresholds.

## Shipping

Ship Monday through Friday. Friday shipments must be labeled for Saturday delivery. All specimens must be labeled with patient's name and collection date. A Viracor Eurofins test requisition form must accompany each specimen. Multiple tests can be run on one specimen. Ship specimens FedEx Priority Overnight® to: Viracor Eurofins, 1001 NW Technology Dr, Lee's Summit, MO 64086.

## Causes for Rejection

Invalid specimen type, inadequate volume, gross hemolysis or gross lipemia, sample not frozen upon receipt.

## Disclaimer

Specimens are approved for testing in New York only when indicated in the Specimen Information field above.

The CPT codes provided are based on Viracor Eurofins' interpretation of the American Medical Association's Current Procedural Terminology (CPT) codes and are provided for general informational purposes only. CPT coding is the sole responsibility of the billing

[https://webdev.viracor-eurofins.com/images/Viracor\\_Eurofins\\_PDFlogo.jpg](https://webdev.viracor-eurofins.com/images/Viracor_Eurofins_PDFlogo.jpg) 1001 NW Technology Drive, Lee's Summit, MO 64086 // (800) 305-5198 // (816) 347-0143 Fax // [info@viracor-eurofins.com](mailto:info@viracor-eurofins.com)

party. Questions regarding coding should be addressed to your local Medicare carrier. Viracor Eurofins assumes no responsibility for billing errors due to reliance on the CPT codes illustrated in this material.

## References

August, KJ, KY Chiang, RM Bostick, WE Flanders, EK Waller, A Langston, D Worthington-White, P Rowland, KF Moore, HJ Khoury and JT Horan. 2011. Biomarkers of immune activation to screen for severe, acute GVHD. *Bone Marrow Transpl* 46: 601-604.

dePablo, R, J Monserrat, E Reyes, D Diaz-Martin, MR Zapata, F Carballo, A de la Hera, A Prieto and M Alvarez-Mon. 2011, Mortality in patients with septic shock correlates with anti-inflammatory but not proinflammatory immunomodulatory molecules. *J Intensive Care Med* 26(2): 125-132.

Sakata, N, M Yasui, T Okamura, M Inoue, K Yumura-Yagi and K Kawa. 2001. Kinetics of plasma cytokines after hematopoietic stem cell transplantation from unrelated donors: the ratio of plasma IL-10/sTNFR level as a potential prognostic marker in severe acute graft-versus-host disease.

Spoetl, T, M Hausmann, F Klebl, A Dirmeier, B Klump, J Hoffmann, H Herfarth, A Timmer, G Rogler. 2007. Serum soluble TNF receptor I and II levels correlate with disease activity in IBD patients. *Inflamm Bowel Dis*. 13(6): 727-732.

Glossop, JR, PT DSawes, NB Nixon, DL Matthey. 2005. Polymorphism in the tumour necrosis factor receptor II gene is associated with circulating levels of soluble tumour necrosis factor receptors in rheumatoid arthritis. *Arthritis Res Ther* 7(6): R1227-1234.

Majewski, S, A Wojas-Pelc, M Malejczyk, E Szymanska and S Jablonska. 1999. Serum levels of soluble TNF alpha receptor type I and the severity of systemic sclerosis. *Acta Derm Venereol* 79(3): 207-210.

Choi SW, P Stiff, K Cooke, JLM Ferrara, T Braun, C Kitko, P Reddy, G Yanik, S Mineishi, S Paczesny, D Hanauer, A Pawarode, E Peres, T Rodriguez, S Smith and JE Levine. 2012. TNF-Inhibition with Etanercept for Graft-versus-Host disease prevention in high risk HCT: Lower TNFRI levels correlate with better outcomes. *Biol Blood Marrow Transplant: Biol Blood Marrow Transpl* 18: 1517-1524.

Paczesny, S, OI Krijanovski, TM Braun, SW Choi, SG Clouthier, R Kuick, DE Misek, KR Cooke, CL Kitko, A Weyand, D Bickley, D Jones, J Whitfield, P Reddy, JE Levine, SM Hanash and JLM Ferrara. 2009. A biomarker panel for acute graft-versus-host disease. *Blood* 113(2): 273-278.

John E. Levine, Brent R. Logan, Juan Wu, Amin M. Alousi, Javier Bolaños-Meade, James L.M. Ferrara, Vincent T. Ho, Daniel J. Weisdorf and Sophie Paczesny (2012) Acute graft-versus-host disease biomarkers measured during therapy can predict treatment outcomes: a Blood and Marrow Transplant Clinical Trials Network study. *Blood*. 119(16):3854-3860.